

# PRO

PROGEAR™

is a line of high-quality **face masks and respirators** made in the United States by Prestige Ameritech.



## PRO INNOVATIONS AND ADVANTAGES

# PRO

PROGEAR™

Manufactured by **PRESTIGE AMERITECH**  
7201 Iron Horse Blvd.  
North Richland Hills, TX 76180

Pro Quality.  
Pro Function.  
Pro Comfort.

# PRO

PROGEAR™

### WAVE-FIT™ NOSEPIECE

All ProGear masks, except the ones with visors, feature our patent-pending Wave-Fit pre-formed nosepiece. This is especially helpful on white masks, as it visually indicates the outside of the masks and prevents users from putting the masks on backwards.



### WINDOWED CARTONS

Another feature unique to ProGear masks is our see-through window on the front of ProGear dispenser boxes. The window lets you know when it is time to replenish.


### ALL PROGEAR MASKS ARE AMERICAN-MADE

Of the leading facemask manufacturers,\* only Prestige Ameritech manufactures masks in the United States.

\*Kimberly-Clark and Cardinal Healthcare do not make masks in the United States

### TAB-LESS EAR-LOOPS

All of Prestige Ameritech's ear-loops are made on state-of-the-art mask machines that leave no waste material hanging from the top or bottom of the mask. Older, less sophisticated mask-manufacturing machines tend to leave these useless appendages.



PRESTIGE AMERITECH'S  
TAB-LESS EAR-LOOP  
MASK



OUR COMPETITION'S  
TABBED EAR-LOOP  
MASK

P: (817) 427-2700 | F: (817) 886-2733  
[www.progearmasks.com](http://www.progearmasks.com)



Made in the  
**USA**

## ANTI-FOG

### CLEAR-VIEW™ ANTI-FOG FACE MASKS

ProGear Clear-View Anti-Fog masks give you three different ways to prevent the fogging of eyeglasses and scopes while wearing a surgical mask. User preference dictates the selection of a tape-strip, foam-strip or film-strip topped mask.



ITEM NO.	DESCRIPTION	QTY
FM80110	Tie-On, Anti-Fog, Tape Nose-Strip	300
FM80120	Tie-On, Anti-Fog, Foam Nose-Strip	300
FM80130	Tie-On, Anti-Fog, Film Nose-Strip	300

ITEM NO.	RATINGS			
	DELTA P	PFE	SPLASH	BFE
FM80110	2.0	98+	No Fluid Resistant Claims	98+
FM80120	2.0	98+	No Fluid Resistant Claims	98+
FM80130	2.0	98+	No Fluid Resistant Claims	98+

## FLUID RESISTANT

### DRY-SHIELD™ PREMIUM FLUID RESISTANT FACE MASKS

Dry-Shield™ Fluid Resistant masks provide the highest fluid resistance of any mask in our line.



In addition to outstanding fluid protection, Dry-Shield masks offer exceptional protection to the wearer with BFE and PFE ratings above 99%.



Dry-Shield masks are available with or without anti-fog visors.



### PRO-80™ FLUID RESISTANT FACE MASKS

Pro-80™ Fluid Resistant masks provide basic fluid resistance against bodily and other fluids at an 80 mm/Hg splash rating.

ITEM NO.	DESCRIPTION	QTY
FM82110	Tie-On, Dry-Shield Premium Fluid Resistant	300
FM82010	Ear-Loop, Dry-Shield Premium Fluid Resistant	400
FM82120	Tie-On, Dry-Shield Premium Fluid Resistant with Anti-Fog Visor	100
FM82020	Ear-Loop, Dry-Shield Premium Fluid Resistant with Anti-Fog Visor	100
FM84110	Tie-On, Pro-80 Standard Fluid Resistant	300
FM84010	Ear-Loop, Pro-80 Standard Fluid Resistant	500

ITEM NO.	RATINGS			
	DELTA P	PFE	SPLASH	BFE
FM82110	3.9	99+	160 mm/Hg	99+
FM82010	3.9	99+	160 mm/Hg	99+
FM82120	3.9	99+	160 mm/Hg	99+
FM82020	3.9	99+	160 mm/Hg	99+
FM84110	2.6	99+	80 mm/Hg	99+
FM84010	2.6	99+	80 mm/Hg	99+

## STANDARD

### PRO-LITE™ SURGICAL MASK

ProGear Pro-Lite Surgical masks are designed to offer superior protection to the wearer while still remaining cool and breathable during long procedures.



ITEM NO.	DESCRIPTION	QTY
FM81110	Tie-On, Surgical - Blue	300

ITEM NO.	RATINGS			
	DELTA P	PFE	SPLASH	BFE
FM81110	2.2	98+	No Fluid Resistant Claims	98+



### PROCEDURE MASK

ProGear's Procedure Masks provide great bacterial protection to the wearer during quick procedures.

ITEM NO.	DESCRIPTION	QTY
FM81020	Ear-Loop, Standard - Blue	500

ITEM NO.	RATINGS			
	DELTA P	PFE	SPLASH	BFE
FM81020	2.2	98+	No Fluid Resistant Claims	98+

## SENSITIVE SKIN

### SKIN-SAVER™ SENSITIVE SKIN MASKS

Skin-Saver™ Sensitive Skin masks are specifically designed for professionals with extra-sensitive skin. Skin-Saver masks contain no colorants, no adhesives, no inks, and no dyes.



ITEM NO.	DESCRIPTION	QTY
FM83110	Tie-On, Sensitive Skin - White	300
FM83010	Ear-Loop, Sensitive Skin - White	500

ITEM NO.	RATINGS			
	DELTA P	PFE	SPLASH	BFE
FM83110	2.0	98+	No Fluid Resistant Claims	98+
FM83010	2.0	98+	No Fluid Resistant Claims	98+

## CUSTOM-PRINTED

### IMPRINTS™ CUSTOM-PRINTED FACE MASK

ProGear Imprints Custom-Printed masks let you get your message out. We can put just about any repeating message, image or graphic on a mask. The minimum order is 100 cases and there is a \$250 one-time plate charge.



ITEM NO.	DESCRIPTION	QTY
FM85120	Tie-On, Customer Selected Print	300
FM85010	Ear-Loop, Customer Selected Print	500
FM85020	Ear-Loop, Isolation, Customer Selected Print	500

ITEM NO.	RATINGS			
	DELTA P	PFE	SPLASH	BFE
FM80110	2.6	99+	No Fluid Resistant Claims	99+
FM80120	2.6	99+	No Fluid Resistant Claims	99+
FM85020	2.6	99+	80 mm/Hg	99+

### PRO-KIDS MASK

ProGear Pro-Kids masks are specially sized for a child's comfort and protection. We use our Imprints™ printing technology to offer many different patterns to encourage children to wear the masks.



ITEM NO.	DESCRIPTION	QTY
FM85030	Ear-Loop, Pediatric	750

ITEM NO.	RATINGS			
	DELTA P	PFE	SPLASH	BFE
FM85030	2.6	99+	No Fluid Resistant Claims	99+

## RESPIRATOR

### N95 PARTICULATE RESPIRATOR

ProGear's N95 Respirator is intended to protect the wearer from airborne particulate materials.



ITEM NO.	DESCRIPTION	QTY
RP88020	N95 Particulate Respirator - White	300

## ISOLATION

### ISOLATION MASKS

ProGear Isolation Masks are designed for use when dealing with patients or materials that are not considered high risk. In addition to standard yellow and green masks, we also offer tissue isolation masks, as well as custom-printed mask designs based on your specifications.



ITEM NO.	DESCRIPTION	QTY
FM81000	Ear-Loop, Tissue Isolation - Yellow	500
FM81010	Ear-Loop, Isolation - Yellow	500
FM81030	Ear-Loop, Isolation - Green	500
FM85020	Ear-Loop, Isolation, Customer Selected Print	500

ITEM NO.	RATINGS			
	DELTA P	PFE	SPLASH	BFE
FM81000	2.0	98+	No Fluid Resistant Claims	98+
FM81010	2.0	98+	No Fluid Resistant Claims	98+
FM81030	2.0	98+	No Fluid Resistant Claims	98+
FM85020	2.6	99+	80 mm/Hg	99+

To find out how to purchase ProGear products, visit our website at [www.progearmasks.com](http://www.progearmasks.com)