



PRO

PROGEAR[®]

RESPIRATORY & PERSONAL PROTECTION PRODUCTS

Manufactured by **PRESTIGE AMERITECH**



- 100% Made in USA**
- WHOLLY DOMESTIC** Raw Materials
- QUALIFYING N95 RESPIRATORS** for Medicare Payment Adjustments*
- CERTIFIED DIVERSE** Native American-Owned Supplier**
- RECOGNIZED** Historically Underutilized Business***
- ESTABLISHED** in 2005
- ISO 13485:2016** Certified
- FDA REGISTERED** Facility
- DESIGNER & BUILDER** of Automated Manufacturing Equipment
- IN-HOUSE** Product Development Team

*Prestige Ameritech's NIOSH-approved ProGear surgical N95 respirators are wholly made in the United States (i.e., based on the Berry Amendment [10 USC § 2533a], the Qualifying N95 Respirators and all of their components are grown, reprocessed, reused, or produced in the United States).

**Certified by the Cherokee Nation Tribal Employment Rights Office (TERO) as an Indian-Owned Company

***Certified by the Texas Comptroller of Public Accounts (CPA) as a State of Texas Historically Underutilized Business (HUB)

PRO[®] PROGEAR[®]

by **PRESTIGE AMERITECH**



RESPIRATORY PROTECTION

Available options include:

- ▶ Earloop, Tie-On, or Horizontal Tie-On Style
- ▶ N95 Respirator and Surgical Masks
- ▶ ASTM Levels 1-3 Fluid Resistance
- ▶ Anti-Fog Foam, Film, or Tape Nosepiece
- ▶ Masks with Eye Shield Visor
- ▶ Breast Cancer Awareness Masks
- ▶ Masks with a See-Through Window
- ▶ Sensitive Skin Masks
- ▶ Pediatric Earloop Masks
- ▶ **Imprints Custom Printed Masks***

**Contact us to learn more about how we can print your logo on any procedure or surgical mask with our in-house digital printing technology.*

MADE IN USA

progearmasks.com

N95 RESPIRATORS



- ▶ Contains only domestic raw materials
- ▶ Qualifies for the Medicare reimbursement for wholly domestic NIOSH-approved surgical N95 respirators
- ▶ ASTM Level 3 fluid barrier protection
- ▶ Cleared by the FDA for use as a surgical mask
- ▶ Filters +98% of particles down to 0.1 microns
- ▶ Large filter area allows for easy breathing and speaking
- ▶ Snugly conforms to many different face shapes
- ▶ Stretch-knit headbands enhance wearer comfort
- ▶ Flat-fold, duckbill style available in standard and comfort

STANDARD SURGICAL N95 RESPIRATORS



RP88010

ProGear N95 Particulate Filter Respirator and Surgical Mask, Small, White (USA)
ASTM Level 3 Protection

300 respirators per case



RP88020

ProGear N95 Particulate Filter Respirator and Surgical Mask, Regular, White (USA)
ASTM Level 3 Protection

300 respirators per case

SURGICAL N95 RESPIRATORS WITH COMFORT STRIP



RP88011

ProGear, Surgical N95 Respirator with Comfort Strip, Small, White (USA)
ASTM Level 3 Protection

300 respirators per case



RP88021

ProGear, Surgical N95 Respirator with Comfort Strip, Regular, White (USA)
ASTM Level 3 Protection

300 respirators per case

All of our masks and respirators are FDA approved and made in America using high quality domestic materials

SURGICAL MASKS



FM80110

Clear-View Anti-Fog Mask, Tape, Tie-On
ASTM Level 1 Protection

300 masks per case



FM82132 Coming Soon!

Dry Shield Fluid Resistant Mask,
Anti-Fog Tape, Tie-On

ASTM Level 3 Protection
300 masks per case



FM80120

Clear-View Anti-Fog Mask, Foam, Tie-On
ASTM Level 1 Protection

300 masks per case



FM81110

Pro-Lite Surgical Mask, Tie-On
ASTM Level 1 Protection

300 masks per case



FM80130

Clear-View Anti-Fog Mask, Film, Tie-On
ASTM Level 1 Protection

300 masks per case



FM83110

Skin-Saver Sensitive Skin Mask, Tie-On
ASTM Level 1 Protection

300 masks per case



FM82110

Dry Shield Fluid Resistant Mask,
Anti-Fog Foam, Tie-On

ASTM Level 3 Protection
300 masks per case



FM83120

Skin-Saver Ultra-Sensitive Skin Mask,
Tie-On

ASTM Level 1 Protection
300 masks per case



FM82120

Dry Shield Fluid Resistant Mask,
Anti-Glare/Anti-Fog Visor, Tie-On

ASTM Level 3 Protection
100 masks per case



FM83130

Skin-Saver Ultra-Sensitive Skin Mask,
Anti-Fog Foam, Tie-On

ASTM Level 1 Protection
300 masks per case



FM82130

Dry Shield Fluid Resistant Mask, Tie-On
ASTM Level 3 Protection

300 masks per case



FM84110

Pro-80 Fluid Resistant Mask, Tie-On
ASTM Level 1 Protection

300 masks per case



FM82131 Coming Soon!

Dry Shield Fluid Resistant Mask,
Anti-Fog Film, Tie-On

ASTM Level 3 Protection
300 masks per case



FM84210

Pro-80 Fluid Resistant Mask,
Horizontal Tie

ASTM Level 1 Protection
300 masks per case

PROCEDURE MASKS

Oxafence® Active Protection Has Arrived

It's time for masks to **inactivate viruses*** rather than just trapping or shielding them.

The newest ProGear ASTM Level 3 Surgical Masks contain Oxafence Active Protection technology. Oxafence inactivates 99.9% of tested viruses* after 5 minutes.

99.9%

INACTIVATION OF
TESTED VIRUSES*



AV82030

Surgical Masks with Oxafence, Earloop, Blue
ASTM Level 3 Protection

500 masks per case

Indications for Use: The ProGear Surgical Mask with Oxafence is a single-use, disposable device used to cover the nose and mouth of the wearer to protect from transfer of microorganisms, body fluids and particulates. The outer layer of the face mask is coated with methylene blue (0.47 micrograms per square cm area).

*The ProGear Surgical Mask with Oxafence inactivates single isolates of the following test viruses after five minutes of contact with the treated surface of the face mask in laboratory (in vitro) tests under a lighting intensity of 500 3 25 lux: Influenza A/PR/8/34/H1N1 (TC adapted, ATCC VR-1469), 3 log or 99.9% inactivation; Betacoronavirus 1 OC43 (ATCC VR-1558), 4 log or 99.99% inactivation; Feline Calicivirus F-9 (ATCC VR- 782), 4 log or 99.99% inactivation. Antiviral efficacy on other viruses has not been evaluated. Correlations between in vitro testing results and any clinical event have not been evaluated.

PROCEDURE MASKS



FM81000

Isolation Mask, Earloop, Yellow

ASTM Level 1 Protection

500 masks per case



FM83010

Skin-Saver Sensitive Skin Mask, Earloop, White

ASTM Level 1 Protection

500 masks per case



FM81001

Metal-Free Mask, Earloop, Tan

Ideal for Behavioral Health and Imaging
ASTM Level 1 Protection

400 masks per case



FM84010

Pro-80 Fluid Resistant Mask, Earloop, Blue

ASTM Level 1 Protection

500 masks per case



FM81020

Procedure Mask, Earloop, Blue

ASTM Level 1 Protection

500 masks per case



FM85012

Mask, Earloop, Pink Ribbon Print, White

ASTM Level 1 Protection

500 masks per case



FM82010

Dry Shield Fluid Resistant Mask, Earloop, Anti-Fog Foam Strip, Blue

ASTM Level 3 Protection

400 masks per case



FM85030

ProKids Pediatric Mask, Earloop, Crayon Print

750 masks per case



FM82020

Dry Shield Fluid Resistant Mask, Earloop, Anti-Glare/ Anti-Fog Visor, Blue

ASTM Level 3 Protection

100 masks per case



FM86010

Clear-View Window Mask, Earloop

ASTM Level 1 Protection

100 masks per case



FM82030

Dry Shield Fluid Resistant Mask, Earloop, Blue

ASTM Level 3 Protection

500 masks per case



FM86011

Clear-View Window Mask, Earloop

ASTM Level 3 Protection

100 masks per case

PRO[®] PROGEAR[®]

by **PRESTIGE AMERITECH**



FACE & EYE PROTECTION

ProGear Face Shields and Disposable Glasses feature lightweight, anti-fog protection from splashing or spraying. All of our shields have optically clear polyester for maximum visibility, and can be worn with or without glasses.

MADE IN USA
progearmasks.com

FACE SHIELDS



FM89010
Face Shield, Full Length,
Foam Top, Elastic
Headband

96 shields per case



FM89030
Face Shield, XL Length,
Foam Top, Elastic
Headband

72 shields per case



FM89011
Face Shield, Full Length,
Loupes Headlight Cut-Out,
Foam Top, Elastic
Headband

96 shields per case



FM89012
Face Shield, Full Length,
11" Drape, Foam Top,
Elastic Headband

96 shields per case

DISPOSABLE GLASSES




FM89110
Pre-Assembled Disposable
Glasses, Assorted Colors

*100 glasses per case
Packed in 10 bags of 10 each*



FM89111
Pre-Assembled Disposable
Glasses, Pink

*100 glasses per case
Packed in 10 bags of 10 each*



FM89112
NEW! Pre-Assembled Disposable Glasses,
ProGear Pack Dispenser

*160 glasses per case
Packed in 4 dispensers of 40 each*

For use with both wall-mounted and flat surface applications.
Dispenser measures 6 in wide by 6.75 in deep by 21 in tall.



PRESTIGE AMERITECH

7201 Iron Horse Blvd.
North Richland Hills, TX 76180
(817) 427-2700

prestigeameritech.com